Hydrofera Blue™

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Descrizione

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Instructions for Use

Description

Hydrofera Bacteriostatic Foam Dressing is a sterile absorptive foam dressing made of Hydrofera Bacteriostatic (PVA) polyvinyl alcohol sponge, Methylene Blue (less than or equal to 0.00025 gr/gr), Crystal Violet (less than or equal to 0.00025 gr/gr).

Indications

Hydrofera Bacteriostatic Foam Wound Dressings are intended as an external dressing for use in local wound management such as pressure ulcers, venous stasis ulcers, arterial ulcers, diabetic ulcers, donor sites, abrasions, lacerations, radiation burns, post-surgical incisions, and other external wounds caused by trauma.

Contraindications

Hydrofera Blue Bacteriostatic Foam Wound Dressings are not indicated for third degree burns.

Application of Hydrofera Blue Dressing

1) Wound Site Preparation

- a. Prior to the initial application of Hydrofera Blue Wound Dressing, the wound should be debrided of any necrotic tissue.
- b. Prior to the application of each Hydrofera Blue Wound Dressing, the wound should be cleansed with normal saline and/or an appropriate cleansing solution.

2) Dressing Application

- a. Open the Hydrofera Blue Foam Wound Dressing package, keeping the dressing on the inside of the package.
- b. Make sure the dressing is larger than the wound area, covering any erythematic tissue. For packing, the dressing may be trimmed to size.

c. Moisten the dressing with sterile saline or sterile water and squeeze out excess.

- d. Place the dressing on or in the wound. The dressing may be folded or layered to ensure the desired fit. (Fill the wound cavity or make sure the dressing is larger than the wound area).
- e. Secure with a secondary dressing to prevent displacement and retain moisture (transparent film or petroleum gauze will suffice).

Dressing Change and Removal

- a. Change the dressing every 1-3 days.
- b. Do not let the dressing completely dry out. (If drying out you can cover it with Vaseline gauze to retain moisture.)
- c. Rehydrate dressing as needed depending on drainage amount.
- d. The dressing should be changed if the dressing turns white.
- e. The dressing should be changed upon strike through drainage.

Precautions and Observations

- 1) This material is not an implant.
- 2) **Infections:** All wounds can be expected to contain some microorganisms; however, if true clinical symptoms should develop, appropriate steps as defined by the attending Physician should be taken to address that infection. Regular evaluation and cleansing of infected wounds should be a common practice.

- 3) **Underlying Causes:** Some wounds such as leg ulcers and pressure ulcers, can be non-healing unless steps are taken to correct the underlying pathology. Hydrofera Blue Bacteriostatic Foam Wound Dressings are designed for local wound management and to be part of the overall management program for these types of wounds.
- 4) **Desiccation:** Hydrofera Blue Bacteriostatic Foam Wound Dressings are designed to maintain a moist wound environment that will facilitate wound healing. If a Hydrofera Blue Dressing is allowed to desiccate, rehydration with sterile saline or sterile water is recommended prior to removal.

Note: Hydrofera Blue Bacteriostatic Foam Wound Dressings have two organic pigments bound to the foam to render it bacteriostatic. The dressing has been tested in vitro and found to be a bacteriostatic barrier for the following bacteria: MRSA, VRE, Staphylococcus, Bacillus subtilis, Yersinia enterocolitica, Staph aureus, Serratia marcescens, E-Coli, Epidermidis. The dressing has not been subjected to a double blind clinical study.*

All safety and effectiveness testing is on file at Hydrofera LLC.

*Studies are currently in progress.

6833 Phillips Industrial Boulevard
Jacksonville, Florida 32256 U.S.A.
Tel: 800-998-8580 • 904-880-1229
Fax: 904-886-9517 • www.invotec.net