Tamponi Nasali

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Indicazioni

Epistassi

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6833 Phillips Industrial Blvd. • Jacksonville, Florida 32256 U.S.A.

Tel: 800-998-8580 / 904-880-1229 • Fax: 904-886-9517 • www.invotec.net

PJ SMITH MEDICAL PRODUCTS MARKETING

18 Yeates Close • Thame OX93AR, UK

○ 0197

Nasal Tampons

DIRECTIONS FOR USE - ENGLISH

Please Read Carefully:

NOTE: THE MEDICAL PROFESSIONAL IS RESPONSIBLE FOR PROPER SURGICAL PROCEDURESANDTECHNIQUES.THEFOLLOWINGINSTRUCTIONSAREFORINFORMATION ONLY. THE ATTENDING SURGEON MUST EVALUATE THE APPROPRIATENESS OF THE PROCEDURE BASED ON HIS/HER OWN MEDICAL TRAINING AND EXPERIENCE.

Indications

The Invotec Nasal Tampons are designed as a nasal packing for controlling nasal hemorrhage. They may also be used for nasal packing following septoplasty, rhinoplasty, or other intra-nasal surgical procedures.

Composition

100% polyvinyl alcohol

Invotec nasal tampon is made from a material with a spongy, porous structure which swells when made wet by a liquid. The nasal tampons are well tolerated by the tissues and have an excellent absorption capacity. The material is not fibrous. In a hydrogenated state it adapts very well to the anatomical state and exerts slight pressure on the neighboring tissues.

How to use it

- Greasing:** In a dry compressed state, half of the nasal tampon (from the tip to halfway down) is greased with an antiseptic cream (PVP iodine salve) over the sides in order to facilitate insertion.
- Insertion:** The nasal tampon is picked up by means of tweezers or by fingers and is carefully pushed into the nose until the whole tampon is inside the nose.
- Expansion:** In this position the tampon is saturated with 0.9% NaCl physiological saline solution by means of an injection needle (+/- 5 ml for the 45 mm tampon and +/- 10 ml for the 80 mm tampon). This shall not be necessary with epistaxis. Depending on the individual situation the tampon shall have to be moistened one or more times per day, in order to prevent desiccation, which would signify some loss of elasticity of the tampon. Do not let nasal tampons remain in position for longer than four days. For extra safety, part of the tampon contains a cord. This cord, secured to the cheek by means of a plaster, is an additional safety factor for the nasal tampons in order to prevent inhalation danger.
- Removal:** Before one removes the nasal tampon, it must be well moistened with 0.9% NaCl physiological saline solution (let it sink in 5 to 10 minutes). On removing the nasal tampons by means of a ventilation tube, it is handy if one first removes the tube, to then carry out the moistening of the tampon through the created hole.

Attention

- The nasal tampons may not be moistened with an alcoholic solution, as then the structure of the material will change.
- The nasal tampons are sterilized.
- The nasal tampons do not contain antibiotics. They do not work bacteriostatically as regards all infections present, nor do they prevent the existence of new infections. However, they shall not strengthen the production of bacteria that cause T.S.S.* to come into being. If an infection already exists, it is recommended that an antimicrobial treatment is carried out before using the nasal tampons. Do not implant this material.
- If the packaging is damaged, the products are no longer suitable for use.
- The nasal tampons can only be used once and may not be re-sterilized.

** Philip M. Tierno, Jr., Ph.D.: Toxic Shock Syndrome Nasal Surgery, Laryngoscope, Volume 97, No. 12, December 1987*

Warranty

Invotec International, Inc. warrants that the product is free from defects in material and workmanship. Invotec will replace or provide a refund for any product found to be defective so long as the product is returned according to the Returned Goods instructions in the Sales Policy. Invotec shall not be liable for any consequential loss, damage or expense directly or indirectly arising from the use of, or inability to use, this product. The foregoing warranty is in lieu of and excludes all other warranties, expressed or implied, however arising, including merchantability, fitness for a particular purpose, against infringement or otherwise. Invotec International neither assumes, nor authorizes any person to assume for it, any other additional liability or responsibility with respect to this product.

— INVOTEC INTERNATIONAL, INC. —

6833 Phillips Industrial Blvd. • Jacksonville, Florida 32256 U.S.A.

Tel: 800-998-8580 / 904-880-1229 • Fax: 904-886-9517 • www.invotec.net

PJ SMITH MEDICAL PRODUCTS MARKETING

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