## Sonda Nasal Para Epistaxis Con Balon Post Stat

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Crown House, Hornbeam Sq. N. • Harrogate HG2 8PB, UK

www.invotec.net

## Post-Stat™ Nasal Balloon Catheter

### Please Read Carefully:

NOTE: THE MEDICAL PROFESSIONAL IS RESPONSIBLE FOR PROPER SURGICAL PROCEDURES AND TECHNIQUES. THE FOLLOWING INSTRUCTIONS ARE FOR INFORMATION ONLY. THE ATTENDING SURGEON MUST EVALUATE THE APPROPRIATENESS OF THE PROCEDURE BASED ON HIS/HER OWN MEDICAL TRAINING AND EXPERIENCE.

### Indications

The Post-Stat nasal balloon catheter was designed as a posterior nasal packing for controlling posterior nasal hemorrhage and as an inter-operative packing to aid in the prevention of fluid aspiration. The Post-Stat catheter allows for both suction and irrigation.

### How to use the Post-Stat Balloon Catheter

- Step 1: Test the cuff before introducing the catheter by inflating with 5 to 10 CC's of sterile saline solution or sterile water.
- Step 2: Lubricate the catheter with sterile water or non-petroleum based lubricant.
- Step 3: If blood clots are present remove the clots from the nose.
- Step 4: Insert the catheter so that the bevelled tip rest in the nasopharynx.
- Step 5: Withdraw the stylet.
- Step 6: Inflate the nasopharyngeal cuff with either sterile water or sterile saline by inserting the syringe in the inflation valve. A 20 CC syringe is provided.  
NOTE: Inflation volumes of 4 to 8 CC's are usually adequate.
- Step 7: Remove the syringe and gently pull the catheter to seat it.
- Step 8: If suction and/or irrigation is to be used set up proper connectors.  
NOTE: Inflation syringe may also be used for introducing irrigation fluids or evacuating blood and irrigation fluids.
- Step 9: Check cuff pressure frequently and relieve cuff pressure when indicated.

### Removal of Catheter

Insert syringe firmly into valve, and evacuate cuff completely. After verification of complete cuff deflation, gently withdraw catheter.

### To Sterilize

If re-sterilization is necessary, STEAM AUTOCLAVING is the recommended method. The Post-Stat catheter should be washed with a mild soap, thoroughly rinsed, wrapped in a lint free wrapper and steam autoclaved at 121° C (250° F) for 30 minutes.

### THIS PRODUCT IS FOR SINGLE USE ONLY.

*Caution: U.S. Federal law restricts this device to sale by or on the order of a physician.*

### Warranty

Invotec International, Inc. warrants that the product is free from defects in material and workmanship. Invotec will replace or provide a refund for any product found to be defective so long as the product is returned according to the Returned Goods instructions in the Sales Policy. Invotec shall not be liable for any consequential loss, damage or expense directly or indirectly arising from the use of, or inability to use, this product. THE FOREGOING WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, HOWEVER ARISING, INCLUDING MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, AGAINST INFRINGEMENT OR OTHERWISE. Invotec International neither assumes, nor authorizes any person to assume for it, any other additional liability or responsibility with respect to this product.

*Post-Stat is a trademark of Invotec International, Inc.*

INVOTEC INTERNATIONAL, INC.

6833 Phillips Industrial Blvd. • Jacksonville, Florida 32256 U.S.A. • Tel: 800-998-8580 / 904-880-1229 • Fax: 904-886-9517

Crown House, Hornbeam Sq. N. • Harrogate HG2 8PB, UK

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Post-Stat est une marque déposée d'Invotec International, Inc.