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**Please Read Carefully:**

NOTE: THE MEDICAL PROFESSIONAL IS RESPONSIBLE FOR PROPER SURGICAL PROCEDURES AND TECHNIQUES. THE FOLLOWING INSTRUCTIONS ARE FOR INFORMATION ONLY. THE ATTENDING SURGEON MUST EVALUATE THE APPROPRIATENESS OF THE PROCEDURE BASED ON HIS/HER OWN MEDICAL TRAINING AND EXPERIENCE.

**Application Instructions**

- Step 1: IMPORTANT...Cleanse and dry skin on nose and application area.
- Step 2: Wipe cleansed area with Isopropyl Alcohol Wipe\* and dry throughtly.
- Step 3: Express edema from application area.
- Step 4: Apply skin preparation\* to application area and **DRY THOROUGHLY...VERY IMPORTANT!**
- Step 5: Apply 1/2 inch wide skin closure tape\* starting at the bottom of the nose being careful to overlap each strip to be assured that all the skin in the application area is covered. **(See figure #1)**
- Step 6: Remove release paper from Dorsum Pad and apply over skin closure tape from Radix to Supratip area if desired. Dorsum Pad adhesive should not contact skin. **(See figure #2)**
- Step 7: Shape the Soft-Form AL splint to fit the application area.
- Step 8: Remove the release paper from the hypoallergenic adhesive on the back of the splint wings and apply foam side to nose over the Dorsum Pad and skin closure tape. Mold the Soft-Form Aluminum Splint to achieve proper fit and pressure.
- Step 9: Splint should be removed in seven to eight days.

\*Included in the complete kit only.

**Removal Instructions**

- Step 1: Remove by cutting wings off with scissors, and gently lifting.
- Step 2: Retaping for three or four days is optional by the physician

**Product Information**

20-10200	Soft-Form AL	Small Size	Single	10/Box
20-10201	Soft-Form AL	Small Size	Kit	5 Kits/Box
20-10210	Soft-Form AL	Regular Size	Single	10/Box
20-10211	Soft-Form AL	Regular Size	Kit	5 Kits/Box
20-10220	Soft-Form AL	Large Size	Single	10/Box
20-10221	Soft-Form AL	Large Size	Kit	5 Kits/Box

**Warranty**

Invotec International, Inc. warrants that the product is free from defects in material and workmanship. Invotec will replace or provide a refund for any product found to be defective so long as the product is returned according to the Returned Goods instructions in the Sales Policy. Invotec shall not be liable for any consequential loss, damage or expense directly or indirectly arising from the use of, or inability to use, this product. **THE FOREGOING WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, HOWEVER ARISING, INCLUDING MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, AGAINST INFRINGEMENT OR OTHERWISE.** Invotec International neither assumes, nor authorizes any person to assume for it, any other additional liability or responsibility with respect to this product.

Soft-Form AL is a trademark of Invotec International, Inc.

**Figure #1** Over-lap skin closure tape **AFTER** skin prep is **DRY**.

**Figure #2** If desired place Dorsum pad on top of skin closure tape **(DO NOT LET DORSUM PAD CONTACT SKIN).**

**Figure #3** Remove release paper from the hypoallergenic adhesive, place molded splint over nose, and adjust for proper fit.

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