Please Read Carefully:

NOTE: THE MEDICAL PROFESSIONAL IS RESPONSIBLE FOR PROPER SURGICAL PROCEDURES AND TECHNIQUES. THE FOLLOWING INSTRUCTIONS ARE FOR INFORMATION ONLY. THE ATTENDING SURGEON MUST EVALUATE THE APPROPRIATENESS OF THE PROCEDURE BASED ON HIS/HER OWN MEDICAL TRAINING AND EXPERIENCE.

Indications

The K-Splint was designed to be used to endoscopic sinus surgery or any other procedure where middle turbinate adhesions or lateralizations are problematic.

How to use the K-Splint

- Step 1: K-splints are made of malleable Teflon®, easy to shape and configure to a patients’ septum and middle turbinate anatomy.
- Step 2: To shape, place the splint on a soft surface; use an instrument such as an elevator to compress the tab area to create a small “cup” to receive the middle turbinate.
- Step 3: Compress the main part of the splint to form angles to allow for variations in the septum.
- Step 4: The K-Splints should be positioned on either side of the septum, with the tab placed lateral to the middle turbinate.
- Step 5: A transeptal suture is used to hold the splints in place, retaining the middle turbinates against the septum.
- Step 6: Indicated sinus procedures may then be performed.
NOTE: Avoid placing the transeptal suture too tightly to help prevent septal perforation.
- Step 7: If the patient complains of postoperative septal pain, remove the splints as soon as possible, otherwise, the splints may be removed at about two weeks postoperatively.

Sterilization

K-splints are supplied sterile for immediate surgical usage.

THIS PRODUCT IS FOR SINGLE USE ONLY.

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician.


Warranty

Invotec International, Inc. warrants that the product is free from defects in material and workmanship. Invotec will replace or provide a refund for any product found to be defective so long as the product is returned according to the Return Goods instructions in the Sales Policy. Invotec shall not be liable for any consequential loss, damage or expense directly or indirectly arising from the use of, or inability to use, this product. THE FOREGOING WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, HOWEVER ARISING, INCLUDING MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, AGAINST INFRINGEMENT OR OTHERWISE. Invotec International neither assumes, nor authorizes any person to assume for it, any other additional liability or responsibility with respect to this product.

<u>Catalog No.</u>	<u>Size</u>
20-10640	Small
20-10642	Medium
20-10644	Large

Symbols and definitions:

	Manufacturer: Indicates the medical device manufacturer, as defined in EU directives 90/385/EEC, 93/42/EEC, and 98/79/EC.
	Authorized representatives in the European Community. This symbol shall be accompanied by the name and address of the authorized representative in the European Community, adjacent to the symbol.
	Use by date: Indicates the date after which the medical device is not to be used.
	Batch code: Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number: Indicates the manufacturer's catalog number so that the medical device can be identified.
	Sterilized using ethylene oxide: Indicates a medical device that has been sterilized using ethylene oxide.
	Do not re-sterilize: Indicates a medical device that is not to be re-sterilized.
	Do not use if package is damaged: Indicates a medical device that should not be used if the package has been damaged or opened.
	Do not re-use: Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Consult instructions for use: Indicates the need for the user to consult instructions for use.
	Caution: Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.

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