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Symbols and definitions:

	Manufacturer: Indicates the medical device manufacturer, as defined in EU directives 90/385/EEC, 93/42/EEC, and 98/79/EC.
	Authorized representatives in the European Community. This symbol shall be accompanied by the name and address of the authorized representative in the European Community, adjacent to the symbol.
	Use by date: Indicates the date after which the medical device is not to be used.
	Batch code: Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number: Indicates the manufacturer's catalog number so that the medical device can be identified.
	Sterilized using ethylene oxide: Indicates a medical device that has been sterilized using ethylene oxide.
	Do not resterilize: Indicates a medical device that is not to be resterilized.
	Do not use if package is damaged: Indicates a medical device that should not be used if the package has been damaged or opened.
	Do not re-use: Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Consult instructions for use: Indicates the need for the user to consult instructions for use.
	Caution: Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	Medical Device Indicated

REV092120

INVOTEC INTERNATIONAL, INC.

6833 Phillips Industrial Blvd. • Jacksonville, Florida 32256 U.S.A.
Tel: 800-998-8580 / 904-880-1229 • Fax: 904-886-9517 • www.invotec.net
MDSS GmbH
Schiffgraben 41 • 30175 Hannover, Germany

DIRECTIONS FOR USE - ENGLISH

Please Read Carefully:

NOTE: THE MEDICAL PROFESSIONAL IS RESPONSIBLE FOR PROPER SURGICALPROCEDURES AND TECHNIQUES. THE FOLLOWING INSTRUCTIONS ARE FOR INFORMATION ONLY. THE ATTENDING SURGEON MUST EVALUATE THE APPROPRIATENESS OF THE PROCEDURE BASED ON HIS/HER OWN MEDICAL TRAINING AND EXPERIENCE.

Indications

The Invotec Nasal Tampons are designed as a nasal packing for controlling nasal hemorrhage. They may also be used for nasal packing following septoplasty, rhinoplasty, or other ntra-nasal surgical procedures.

Composition

100% polyvinyl alcohol

Invotec nasal tampon is made from a material with a spongy, porous structure which swells when made wet by a liquid. The nasal tampons are well tolerated by the tissues and have an excellentabsorption capacity. The material is not fibrous. In a hydrogenated state it adapts very well to the anatomical state and exerts slight pressure on the neighboring tissues.

How to use it

- Greasing: In a dry compressed state, half of the nasal tampon (from the tip to halfway down) is greased with an antiseptic cream (PVP iodine salve) over the sides in order to facilitate insertion.
- Insertion: The nasal tampon is picked up by means of tweezers or by fingers and is carefully pushed into the nose until the whole tampon is inside the nose.
- Expansion: In this position the tampon is saturated with 0.9% NaCl physiological saline solution by means of an injection needle (+/- 5 ml for the 45 mm tampon and +/- 10 ml for the 80 mm tampon). This shall not be necessary with epistaxis. Depending on the individual situation the tampon shall have to be moistened one or more times per day, in order to prevent desiccation, which would signify some loss of elasticity of the tampon. Do not let nasal tampons remain in position for longer than four days. For extra safety, part of the tampon contains a cord. This cord, secured to the cheek by means of a plaster, is an additional safety factor for the nasal tampons in order to prevent inhalation danger.
- Removal: Before one removes the nasal tampon, it must be well moistened with 0.9% NaCl physiological saline solution (let it sink in 5 to 10 minutes). On removing the nasal tampons by means of a ventilation tube, it is handy if one first removes the tube, to then carry out the moistening of the tampon through the created hole.

Attention

- The nasal tampons may not be moistened with an alcoholic solution, as then the structure of the material will change.
 - The nasal tampons are sterilized.
 - The nasal tampons do not contain antibiotics. They do not work bacteriostatically as regards all infections present, nor do they prevent the existence of new infections. However, they shall not strengthen the production of bacteria that cause T.S.S.* to come into being. If an infection already exists, it is recommended that an antimicrobial treatment is carried out before using the nasal tampons. Do not implant this material.
 - If the packaging is damaged, the products are no longer suitable for use.
 - The nasal tampons are intended for SINGLE USE ONLY. Reuse of this device may expose the patient to infection or contamination risks. Once the seal of the sterile packaging has been torn open its contents will not be taken back by the manufacturer.
- * Philip M. Tierno, Jr., Ph.D.: *Toxic Shock Syndrome Nasal Surgery, Laryngoscope, Volume 97, No.12, December 1987*

Warranty

Invotec International, Inc. warrants that the product is free from defects in material and workmanship. Invotec will replace or provide a refund for any product found to be defective so long as the product is returned according to the Returned Goods instructions in the Sales Policy. Invotec shall not be liable for any consequential loss, damage or expense directly or indirectly arising from the use of, or inability to use, this product. The foregoing warranty is in lieu of and excludes all other warranties, expressed or implied, however arising, including merchantability, fitness for a particular purpose, against infringement or otherwise. Invotec International neither assumes, nor authorizes any person to assume for it, any other additional liability or responsibility with respect to this product.

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