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**Symbols and definitions:**

	Manufacturer: Indicates the medical device manufacturer, as defined in EU directives 90/385/EEC, 93/42/EEC, and 98/79/EC.
	Authorized representatives in the European Community. This symbol shall be accompanied by the name and address of the authorized representative in the European Community, adjacent to the symbol.
	Use by date: Indicates the date after which the medical device is not to be used.
	Batch code: Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number: Indicates the manufacturer's catalog number so that the medical device can be identified.
	Sterilized using ethylene oxide: Indicates a medical device that has been sterilized using ethylene oxide.
	Do not resterilize: Indicates a medical device that is not to be resterilized.
	Do not use if package is damaged: Indicates a medical device that should not be used if the package has been damaged or opened.
	Do not re-use: Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Consult instructions for use: Indicates the need for the user to consult instructions for use.
	Caution: Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.

**Please Read Carefully:**

NOTE: THE MEDICAL PROFESSIONAL IS RESPONSIBLE FOR PROPER SURGICAL PROCEDURES AND TECHNIQUES. THE FOLLOWING INSTRUCTIONS ARE FOR INFORMATION ONLY. THE ATTENDING SURGEON MUST EVALUATE THE APPROPRIATENESS OF THE PROCEDURE BASED ON HIS/HER OWN MEDICAL TRAINING AND EXPERIENCE.

**Indications**

When a patient is a poor candidate for surgical correction of a septal perforation, and the nature of the perforation, that is size and location, and septal anatomy are appropriate, the attending surgeon may consider closure with a septal button. Generally, anterior septal perforations that are 2cm. or less in maximum diameter are suitable for closure with a septal button. This product is for single use only. Reuse of this device may expose the patient to infection or contamination risks. Once the seal of the sterile packaging has been torn open, its contents will not be taken back by the manufacture.

**Contraindications**

Septal deformities are relative contraindications for septal button placement, as the flange on the convex side of the deformity may increase nasal obstruction. Relative contraindications may also include active local nasal infection as the septal button may prevent adequate treatment.

**How to use the Septal Button**

Step 1: After adequate topicalization of the nose, 1% lidocaine with 1:100,000 epinephrine is infiltrated about the margins of the perforation. Light sedation may be required in some cases.

Step 2: Custom trim the flanges of the button to permit 3-4mm. of overlap with respect to the margins of the perforation.

Step 3: Choose one of the following two methods of insertion:

**METHOD 1:**

Collapse one flange and introduce the button into the perforation, then gently pull it back through from the opposite end.

**METHOD 2:**

Sew the edges of one flange together in a rosette manner with a suture, tie suture and leave long. Pass the long end of the suture through the perforation and gently guide the button in place using a hemostat from the back side to apply gentle pressure. Cut and remove the suture allowing the sewn flange to open flat.

**NOTE: CARE MUST BE TAKEN SO THAT THE FLANGES OF THE BUTTON DO NOT ABUTT THE JUNCTION OF THE SEPTUM WITH THE UPPER LATERAL CARTILAGE, AS THIS IS A PARTICULARLY ANNOYING SENSATION FOR THE PATIENT.**

**Removal**

Device may be used up to twenty-nine (29) days. Remove button by carefully cutting the center post with a scalpel or microscissors and removing each piece of the button through separate nostrils.

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician.

**Warranty**

Invotec International, Inc. warrants that the product is free from defects in material and workmanship. Invotec will replace or provide a refund for any product found to be defective so long as the product is returned according to the Returned Goods instructions in the Sales Policy. Invotec shall not be liable for any consequential loss, damage or expense directly or indirectly arising from the use of, or inability to use, this product. THE FOREGOING WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, HOWEVER ARISING, INCLUDING MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, AGAINST INFRINGEMENT OR OTHERWISE.

Invotec International neither assumes, nor authorizes any person to assume for it, any other additional liability or responsibility with respect to this product.

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