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Symbols and definitions:

	Manufacturer: Indicates the medical device manufacturer, as defined in EU directives 90/385/EEC, 93/42/EEC, and 98/79/EC.
	Authorized representatives in the European Community. This symbol shall be accompanied by the name and address of the authorized representative in the European Community, adjacent to the symbol.
	Use by date: Indicates the date after which the medical device is not to be used.
	Batch code: Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number: Indicates the manufacturer's catalog number so that the medical device can be identified.
	Do not use if package is damaged: Indicates a medical device that should not be used if the package has been damaged or opened.
	Do not re-use: Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Consult instructions for use: Indicates the need for the user to consult instructions for use.
	Caution: Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.

Please Read Carefully:

NOTE: THE MEDICAL PROFESSIONAL IS RESPONSIBLE FOR PROPER SURGICAL PROCEDURES AND TECHNIQUES. THE FOLLOWING INSTRUCTIONS ARE FOR INFORMATION ONLY. THE ATTENDING SURGEON MUST EVALUATE THE APPROPRIATENESS OF THE PROCEDURE BASED ON HIS/HER OWN MEDICAL TRAINING AND EXPERIENCE.

**Application Instructions**

Step 1: IMPORTANT! Cleanse and dry skin on nose and application area.

Step 2: Wipe cleansed area with Isopropyl Alcohol Wipe\* and dry thoroughly.

Step 3: Express edema from application area.

Step 4: Apply skin preparation\* to application area and DRY THOROUGHLY (Extremely IMPORTANT!)

Step 5: Apply 1/2 inch wide skin closure tape\* starting at the bottom of the nose being careful to overlap each strip to be assured that all the skin in the application area is covered. (See Figure #1)

Step 6: Remove release paper from Dorsum Pad and apply over skin closure tape from Radix to Supratip area if desired. Dorsum Pad adhesive should not contact skin. (See Figure #2)

Step 7: Shape the Soft-Form AL splint to fit the application area.  
Step 8: Remove the release paper from the hypoallergenic adhesive on the back of the splint wings and apply foam side to nose over the Dorsum Pad and skin closure tape. Mold the Soft-Form Aluminum Splint to achieve proper fit and pressure. (See Figure #3)

Step 9: Splint should be removed in seven to eight days.

\*Included in the complete kit only.

**Removal Instructions**

Step 1: Remove by cutting wings off with scissors, and gently lifting.

Step 2: Retaping for three or four days is optional by the physician

**Warranty**

Invotec International, Inc. warrants that the product is free from defects in material and workmanship. Invotec will replace or provide a refund for any product found to be defective so long as the product is returned according to the Returned Goods instructions in the Sales Policy. Invotec shall not be liable for any consequential loss, damage or expense directly or indirectly arising from the use of, or inability to use, this product. THE FOREGOING WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, HOWEVER ARISING, INCLUDING MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, AGAINST INFRINGEMENT OR OTHERWISE. Invotec International neither assumes, nor authorizes any person to assume for it, any other additional liability or responsibility with respect to this product.

Soft-Form AL™ is a trademark of Invotec International, Inc.

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Yumuşak-Yapılı AL bir Invotec International, Inc. ticari markasıdır.