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Symbols and definitions:

	Manufacturer: Indicates the medical device manufacturer, as defined in EU directives 90/385/EEC, 93/42/EEC, and 98/79/EC.
	Authorized representatives in the European Community. This symbol shall be accompanied by the name and address of the authorized representative in the European Community, adjacent to the symbol.
	Use by date: Indicates the date after which the medical device is not to be used.
	Batch code: Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number: Indicates the manufacturer's catalog number so that the medical device can be identified.
	Sterilized using ethylene oxide: Indicates a medical device that has been sterilized using ethylene oxide.
	Do not re-sterilize: Indicates a medical device that is not to be re-sterilized.
	Do not use if package is damaged: Indicates a medical device that should not be used if the package has been damaged or opened.
	Do not re-use: Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Consult instructions for use: Indicates the need for the user to consult instructions for use.
	Caution: Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	Medical Device Indicated

REV 092120

Please Read Carefully:

NOTE: THE MEDICAL PROFESSIONAL IS RESPONSIBLE FOR PROPER SURGICAL PROCEDURES AND TECHNIQUES. THE FOLLOWING INSTRUCTIONS ARE FOR INFORMATION ONLY. THE ATTENDING SURGEON MUST EVALUATE THE APPROPRIATENESS OF THE PROCEDURE BASED ON HIS/HER OWN MEDICAL TRAINING AND EXPERIENCE.

Indications

The Ultra-Stat Catheter was designed as a nasal packing for controlling nasal hemorrhage. Ultra-Stat may also be used for nasal packing following septoplasty, rhinoplasty, or other intra-nasal surgical procedures.

How to use the Ultra-Stat Balloon Catheter

Step 1: Test inflate the cuff before introducing the catheter. Each inflation valve is marked to indicate its corresponding cuff by volume.

Step 2: Lubricate the catheter with sterile water or non-petroleum based lubricant.

Step 3: Remove blood clots from nose.

Step 4: Insert the catheter so that the beveled tip rest in the nasopharynx.

Step 5: Inflate the nasopharyngeal cuff with either sterile water or sterile saline by inserting the syringe in the valve marked 10 CC. A 20 CC syringe is provided.

Step 6: Remove the syringe and gently pull the catheter to seat it.

Step 7: Inflate the anterior balloon in the same manner using up to 30 CC's of sterile solution.

NOTE: Inflation of 10 to 25 CC's are usually adequate.

Step 8: Check both cuff pressures frequently and relieve cuff pressures when indicated.

NOTE: A small, well lubricated suction catheter may be inserted through the airway to remove secretions and/or blood clots.

Removal of Catheter

Insert syringe firmly into each valve, in turn, and evacuate each cuff completely. After verification of complete cuff deflation, gently withdraw catheter.

General Handling

Any contamination of the product must be avoided.

The product must be kept in its sealed protective packaging.

Examine packaging for damage before opening since damaged packaging may impair the sterility of the product.

Do not open the protective until just before implantation of the product. In addition, the product must be visually checked for damage.

This product is intended for single use only. Reuse of this device may expose patient to infection or contamination risks. Once the seal of the sterile packaging has been torn open its contents will not be taken back by the manufacturer.

THIS PRODUCT IS FOR SINGLE USE ONLY.

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician.

Warranty

Invotec International, Inc. warrants that the product is free from defects in material and workmanship. Invotec will replace or provide a refund for any product found to be defective so long as the product is returned according to the Returned Goods instructions in the Sales Policy. Invotec shall not be liable for any consequential loss, damage or expense directly or indirectly arising from the use of, or inability to use, this product. THE FOREGOING WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, HOWEVER ARISING, INCLUDING MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, AGAINST INFRINGEMENT OR OTHERWISE. Invotec International neither assumes, nor authorizes any person to assume for it, any other additional liability or responsibility with respect to this product.

Ultra-Stat is a trademark of Invotec International, Inc.

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Ultra-Stat Bir Invotec International, Inc. markasıdır.