Leer atentamente:

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Please Read Carefully:

NOTE: THE MEDICAL PROFESSIONAL IS RESPONSIBLE FOR PROPER SURGICAL PROCEDURES AND TECHNIQUES. THE FOLLOWING INSTRUCTIONS ARE FOR INFORMATION ONLY. THE ATTENDING SURGEON MUST EVALUATE THE APPROPRIATENESS OF THE PROCEDURE BASED ON HIS/HER OWN MEDICAL TRAINING AND EXPERIENCE.

THIS PRODUCT IS INTENDED FOR SINGLE USE ONLY.

Reuse of this device may expose the patient to infection or contamination risks.

Application Instructions

Step 1: IMPORTANT! Cleanse and dry skin on nose and application area.

Step 2: Wipe cleansed area with Isopropyl Alcohol Wipe* and dry thoroughly.

Step 3: Express edema from application area.

Step 4: Apply skin preparation* to application area and DRY THOROUGHLY (Extremely IMPORTANT!)

Step 5: Apply 1/2 inch wide skin closure tape* starting at the bottom of the nose being careful to overlap each strip to be assured that all the skin in the application area is covered. (See Figure #1)

Step 6: Remove release paper from Dorsum Pad and apply over skin closure tape from Radix to Supratip area if desired. Dorsum Pad adhesive should not contact skin. (See Figure #2)

Step 7: Shape the Soft-Form™ AL splint to fit the application area.

Step 8: Remove the release paper from the hypoallergenic adhesive on the back of the splint wings and apply foam side to nose over the Dorsum Pad and skin closure tape. Mold the SoftForm™ Aluminum Splint to achieve proper fit and pressure. (See Figure #3)

Step 9: Splint should be removed in seven to eight days. *Included in the complete kit only.

Removal Instructions

Step 1: Remove by cutting wings off with scissors, and gently lifting.

Step 2: Retaping for three or four days is optional by the physician

Disposal

When disposing the product, take all steps necessary to avoid risk of injury and infection. Contaminated products must be disposed of as hazardous waste and handled so as to avoid contamination of third parties.

Reporting

Any serious incident that has occurred in relation to the device should be reported to Invotec International, Inc, MDSS GmbH, and thew competent authority of the Member State in which the user and/or patient is established.

Warranty

Invotec International, Inc. warrants that the prAoduct is free from defects in material and workmanship. Invotec will replace or provide a refund for any product found to be defective so long as the product is returned according to the Returned Goods instructions in the Sales Policy. Invotec shall not be liable for any consequential loss, damage or expense directly or indirectly arising from the use of, or inability to use, this product. THE FOREGOING WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, HOWEVER ARISING, INCLUDING MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, AGAINST INFRINGEMENT OR OTHERWISE. Invotec International neither assumes, nor authorizes any person to assume for it, any other additional liability or responsibility with respect to this product.

Soft-Form™ AL is a trademarks of Invotec International, Inc.

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Soft-Form™ AL, Invotec International, Inc.'in ticari bir markasıdır.