Leer atentamente:

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Table with 2 columns: Symbol and Description. Symbols include factory, EC REP, hourglass, LOT, REF, STERILE EO, no re-sterilize, damaged package, no reuse, no reuse symbol, open book with 'i', and warning triangle. Descriptions are in multiple languages (ES, PT, FR, IT, DE, NL, TK) and explain the meaning of each symbol.

REV 061021

INVOTEC INTERNATIONAL, INC.

6833 Phillips Industrial Blvd. • Jacksonville, Florida 32256 U.S.A. Tel: 800-998-8580 / 904-880-1229 • Fax: 904-886-9517 • www.invotec.net

MDSS GmbH

Schiffgraben 41 • 30175 Hannover, Germany

Please Read Carefully:

NOTE: THE MEDICAL PROFESSIONAL IS RESPONSIBLE FOR PROPER SURGICAL PROCEDURES AND TECHNIQUES. THE FOLLOWING INSTRUCTIONS ARE FOR INFORMATION ONLY. THE ATTENDING SURGEON MUST EVALUATE THE APPROPRIATENESS OF THE PROCEDURE BASED ON HIS/HER OWN MEDICAL TRAINING AND EXPERIENCE.

Indications

1. COAT THE PACK: Coat the leading edge of the ear pack with an antibiotic ointment or lubricant for ease of insertion. If applicable, trim drawingstring to desired length. Ear Pack may be trimmed to size if required.

2. INSERT: Grasp the ear pack with forceps. Hold the external canal open with a speculum and advance the pack into the external canal. Gelfoam or antibiotic ointment may be inserted first. Instill antibiotic ointment generously around the side of the pack, leaving the proximal end uncoated for hydration.

3. EXPANSION: Hydrate the outer end of the ear pack with saline, sterile water or antibiotic drops until fully expanded.

4. REMOVAL: Re-hydrate the pack until saturated (surgeons recommend allowing the pack to remain in place for 5 - 10 minutes after re-hydrate to facilitate removal). Use forceps to grasp the pack and gently withdraw.

- Surgeons recommend instructing the patient to instill antibiotic drops onto the pack at least twice a day.

Important

PVA Ear Packing is Single Use Only. Reuse of this device may expose patient to infection or contamination risks. once the seal of the sterile packaging has been torn open its contents will not be taken back by the manufacturer.

BE SURE TO RE-HYDRATE THE PACK BEFORE REMOVAL.

Disposal: When disposing the product, take all steps necessary to avoid risk of injury or infection. Contaminated products must be disposed of as hazardous waste and handled to avoid contamination of third parties.

Note

Any serious incident that has occurred in the relation to the device should be reported to Invotec international, Inc. MDSS GmbH, and the competent authority of the Member State in which the user and/or patient is established.

Warranty

Invotec International, Inc. warrants that the product is free from defects in material and workmanship. Invotec will replace or provide a refund for any product found to be defective so long as the product is returned according to the Returned Goods instructions in the Sales Policy. Invotec shall not be liable for any consequential loss, damage or expense directly or indirectly arising from the use of, or inability to use, this product. THE FOREGOING WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, HOWEVER ARISING, INCLUDING MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, AGAINST INFRINGEMENT OR OTHERWISE. Invotec International neither assumes, nor authorizes any person to assume for it, any other additional liability or responsibility with respect to this product.

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