INVOTEC®

VENTILATION TUBES

General handling

Any contamination of the product must be avoided. The product must be kept in its sealed protective packaging. Examine packaging for damage before opening since damaged packaging may impair the sterility of the product. Do not open the protective packaging until just before implantation of the product. In addition, the product must be visually checked for damage. This product is intended for single us only. Once the seal of the sterile packaging has been torn open its contents will not be taken back by the manufacturer.

Packaging and sterility

The packaging of the product consists of Sterile packaging (primary packaging) The packaging meets the provisions of the European Standards. Intact packaging protects the product form environmental influences and ensures sterile storage.

Handling of sterile packaging

Please ensure that the relevant aseptic instructions are complied with when removing the product from the packaging.

Re-sterilization

Re-sterilization after expiry of the use-by date as well as re-sterilization of unsterile uncontaminated products is prohibited. Manufacturer and distributor will assume no liability for products that have been re-sterilized by the user.

Storage

Store in a dry environment at room temperature. After expiry of the use-by date the product may no longer be used.

Product description

Ventilation tube for ventilation of the middle ear. The tube should be implanted via a small incision in the eardrum (myringotomy). The tube has a central lumen of varying diameter. The actual size to used should be decided by the surgeon based on the respective indication. The tube facilitates the removal of secretions (middle ear effusion) and relieving the auditory tube (Eustachian tube).

Indications

- Tympanum (Middle Ear) effusion (secretory otitis media, glue ear)
- Chronic Eustachian tube related middle ear ventilation problems

Contraindications

Allergic reactions to the material used.

Side-Effects

- In rare cases otitis media is possible after myringotomy; in such cases the patient can be treated with either a topical or systemic antibiotic.
- Insertion of an otologic ventilation tube very rarely results in a potentially permanent eardrum perforation in which in turn can be closed through tympanoplastic surgery, if necessary.
- The tube may also fall into the middle ear. In this case the tube can be surgically removed.
- Usually the tube will migrate spontaneously after a prolonged period, when it simply falls out onto the outer ear. If this does not occur then the tube can be surgically removed, though this is not normally required.
- If strong viscid secretions are produced (secretory otitis media) the tube may become clogged. If it is not possible to free the tube from these clogging secretions then the clogged tube can be removed and replaced by a new tube.
- Premature migration of the ventilation tube may occur. If this occurs, new otologic ventilation tube may have to be inserted after the reasons for the early migration have been identified.
- There is the possibility of hyalinization and subsequent calcification (myringosclerosis / tympanosclerosis of the ear drum.
- Improper insertion may cause subluxation of the ossicular chain.
- The existence of anatomical anomalies may lead to injury of other middle ear structures including the facial nerve, tensor tympani or stapedius muscles, or temporal bone related vital structures such as the middle ear, carotid artery or jugular vein.

Precautions

While the tube remains inserted in the eardrum the patient need be advised that water or contaminated fluid can flow through the tube into the middle ear, since this may result in a secondary infection. This can occur as a result of bathing, washing the hair, swimming, etc. To avoid such infection the use of tight fitting ear plugs is recommended. The patient needs to be aware that significant pressure variations while the tube remains in situ, e.g. from deep water diving, may result in injured ear and possibly hearing injury. All ventilation tubes from Invotec International are MRI compatible to 4.5T.

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