

Hemostat Kit

Porcine Gelatin Hemostat Powder

The Invotec Porcine Hemostat Kit contains all the key components to prepare and deliver a hemostat porcine gel. Thrombin may be added to enhance the hemostatic action of the porcine-based gel or putty.



Hemostat Kit Includes

- 2 each 10cc Luer Loc syringes
- 1 16g, 1½", Hypodermic transfer needle
- 17g, 4" (15cm) flexible cannula, Luer Loc hub
- 1 gram Gelitan Powder in sterile mixing tray

Gelitan Packing & Kits

Injectable Crosslinked Hyaluronic Acid Nasal Packing and Stent

Gelita-Spon™ Gel Hemostat



Catalog No.	Description	Qty/Box
EQP-01K	Hemostat Sterile Wrapped Kit, includes: 2 each 10cc Luer Loc syringes 1 16g, 1½", Hypodermic transfer needle 17g, 4" (15cm) flexible cannula, Luer Loc hub 1 gram Gelitan Powder in sterile mixing tray	5


Gelatin Powder



Catalog No.	Description		Qty/Box
EQP-01	Equi-Spon® Powder, 1 gram	Sterile	4

PR-00005C	17g, 4" flexible instillation cannula, Luer Loc hub	Sterile	1
-----------	---	---------	---

Gelatin Sponges and Film USP



Catalog No.	Description	Qty/Box
EQU60-20-07	Sponge, 60 x 20 x 7 mm, Size 12-7	10
EQU125-80-10	Sponge, 125 x 80 x 10 mm, Size 100	5
EQU80-50-10	Sponge, 80 x 50 x 10 mm, Size 50	10
EQU80-30-00	Cylindrical Tampon, 80 x 30 mm	5
EQU200-70-05	Film 200 x 70 x 5.0 mm	10

 **INVOTEC
INTERNATIONAL, INC.®**

6833 PHILLIPS INDUSTRIAL BLVD.
JACKSONVILLE, FL 32256

REV 0817

800 998 8580 • www.invotec.net