

Micro-Needle

Precision, Exquisite Electrosurgery



The Invotec™ Micro-Needle is designed for clean, precise electrosurgical dissection during fine aesthetics procedures. The insulated tungsten needle delivers the generated wave form to the exquisitely small needle dissector tip while retarding eschar and maintaining the finest precision tip.

Using extremely low wattage, driving reduced opportunity for tissue necrosis, the Invotec Micro-Needle provides precision cutting and cautery with minimal plume. during many surgical procedures including:

- | | |
|-------------------------------|---------------------------------------|
| Blepharoplasty | Cleft Lip Revision |
| Rhinoplasty | Pediatric Plastics |
| Neck dissection | Bicoronal Flap |
| Otoplasty | Hand Surgery |
| Tram Flap | Facial Plastic Procedure |
| Incising Temporalis
Fascia | Forehead Corrugator
Muscle Release |

Ordering Information

Catalog No.	Description	Quantity
29-10000	Micro-Needle, Straight, 5.5cm (2.2"), blue, insulated, exposed Tungsten tip 3mm, sterile	1